



 **FRAXFACE**  
FRACTIONAL RADIOFREQUENCY

ECLERIS   
INNOVATION VALUE



3 Programs  
3 Energy Levels  
Deep stimulation effects


3 Effects:  
• ablation (micro ablation)  
• coagulation  
• heating

## FRACTIONAL RF

FraxFace uses radiofrequency energy through an applicator with 64 micro-electrode points of contact with the skin. Each point in contact with the epidermis generates different levels of thermal stimulation at different skin depths, with the fundamental characteristic of leaving the surrounding tissue unaffected. These areas will be the ones responsible for the quick recovery of the skin.

The energy applied has a pyramid pattern, with minimal disruption in the epidermis and wider effect in the deeper layers, so that thermal effect is localized at the dermis-epidermis union. The different programs allow different layers to be reached with full control.

Recovery, renovation and rejuvenation are produced from the deeper layers with minimal and temporal detriment to epidermis.



■ ABLATION ■ COAGULATION ■ HEATING

### Indications:

- Elimination of fine lines
- Scars and pores
- Acne scars
- Skin irregularities
- Photoaging
- Lack of tone and radiance
- Hypertrophic scars
- Stretch marks
- Keratosis

Rejuvenation on face, neck, cleavage and hands

# FRAXFACE



**ECLERIS**   
INNOVATION VALUE  
[www.ecleris.com](http://www.ecleris.com)  
[info@ecleris.com](mailto:info@ecleris.com)

Ecleris reserves the right to modify this product without previous notice.