

EVERLUX MINI

SUPER BRIGHT LED LIGHT SOURCE

LED 80W 50.000, HOURS LED LIFE



BRIGHT, WHITE AND POWERFUL

EVERLUX MINI

SUPER BRIGHT LED LIGHT SOURCE



Interchangeable fiber optic connection.
Compatible with main brands



TECHNICAL SPECIFICATIONS

SERIES	LED LIGHT SOURCE
LED	80 W
LIGHT COLOR TEMPERATURE	5900° K
POWER SUPPLY VOLTAGE	90 - 264 V
N° OF PHASES	1 (SINGLE-PHASE)
FREQUENCY	50 / 60Hz
POWER CONSUMPTION	175 WB
CLASS	I
TYPE	B
LED SYSTEM SERVICE LIFE	50.000 HOURS APPROX.
REFRIGERATION	FORCED THROUGH A FAN
OVER TEMP. PROTECTION	AUTOMATIC LED SYSTEM POWER-OFF
DIMMER	ELECTRONIC
FIBER OPTIC CONNECTOR	INTERCHANGEABLE, COMPATIBLE WITH MAIN BRANDS
DIMENSIONS	202 X 290 X 95mm (width x depth x height)
WEIGHT	2.5 Kg.



www.ecleris.com

ECLERIS
INNOVATION VALUE