


EVERLUX PRO

SUPER BRIGHT LED LIGHT SOURCE

LED 80W 50.000, HOURS LED LIFE



Integrated connector
for optional camera

BRIGHT, WHITE AND POWERFUL



EVERLUX PRO

SUPER BRIGHT LED LIGHT SOURCE



Turret for fiber optic connection compatible with main brands.



TECHNICAL SPECIFICATIONS

SERIES	LED LIGHT SOURCE
LED	80 W
LIGHT COLOR TEMPERATURE	5900°K
POWER SUPPLY VOLTAGEN°	90 - 264 V
OF PHASES	1 (SINGLE-PHASE)
FREQUENCY	50 / 60Hz
POWER CONSUMPTION	175 W
CLASS	I
TYPE	B
LED SYSTEM SERVICE LIFE	50.000 HOURS APPROX.
REFRIGERATION	FORCED THROUGH A FAN
OVER TEMP. PROTECTION	AUTOMATIC LED SYSTEM POWER-OFF
DIMMER	ELECTRONIC
FIBER OPTIC CONNECTOR	BUILT-IN REVOLVING TURRET WITH MULTIPLE ADAPTERS
DIMENSIONS	202 X 290 X 95mm (width x depth x height)
WEIGHT	3.0 Kg.
SERIES	VIDEO CAMERA
RESOLUTION	PAL 625 LINES / NTSC 525 LINES
CCD	1/3"
SHUTTER SPEED	UP TO 1/100.000 SEC.
MAX. ILLUMINATION	2 LUX F1.4
VIDEO OUTPUT	COMPOSITE VIDEO (CVBS) (1 Vpp/75 Ω), BNC S-VIDEO (Y-C), JACK MINI-DIN MONO JACK 3.5mm
ENDODIGI CONNECTOR	WHITE BALANCE
ALL CONTROLS	ANTI-MOIRÉ (FOR FLEXIBLE SCOPES)
ON HANDPIECE	WINDOW PHOTO



www.ecleris.com

ECLERIS
INNOVATION VALUE